



# Leak Test Apparatus CL 707 & CL 708



## • Introduction

**CHEMILINE Leak Test Apparatus CL 707 & CL 708** is very sturdy in operation ideal for checking leakage of parked strip, blisters & small sachets. Meets USP and package practice specifications. Programmable vacuum level & hold time. Slow venting for better result. Validation Print out It is available with 150 mm dia jar, 200 mm dia jar and 300 mm dia jar.

CL 707



CL 708



## • Specifications

### CL 707

- Ideal for checking leakage of parked strips blisters and small sachets.
- Meets USP & package practice specifications.
- Programmable vacuum level & hold time.
- Slow venting for better result.

### CL 708

- Fully automatic and programmable timer and vacuum (Up to -560mm of Hg)
- Automatic release of vacuum at the end of the test
- Automatic control of vacuum level, by oil free vacuum pump
- Polycarbonate vacuum desiccator with instrument
- GMP model – STAINLESS STEEL body
- Compact size makes the instrument portable and cute
- Slow release of vacuum at the end of the test
- Clear 4 x 20 LCD display

**Labline Technologies Pvt. Ltd.**